



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Antithrombin

Concentration: 2.42 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.413 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M sodium
citrate/ 0.15 M NaCl/ pH 8.3

Storage: < -60 °C

Expiration Date: 05/2006

PO: N/A

The Bovine Antithrombin is >95% pure on SDS-PAGE gels and shows no reduction upon incubation with 2-mercaptoethanol. All bovine material is purified from U.S. cows < 30 months of age.

Thank You,
Quality Assurance

Sample