



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Bovine Fibrinogen, Plasminogen Depleted

**Concentration:** 46.72 mg/ml

**Aliquots:** 1 x 1.00 g

**Volume/aliquot:** 21.404 ml

**Activity:** 100.00 % Clottable

**Total protein:** 1.00 g

**Buffer:** 20 mM sodium citrate-HCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 05/2009

**PO:** N/A

The Bovine Fibrinogen, Plasminogen depleted is >95% pure by SDS-PAGE and is greater than 95% clottable. The plasminogen was depleted using a lysine-sepharose column. Protein must be thawed in a 37oC water bath, aliquot and freeze the unused portion at <-60oC. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,  
**Quality Assurance**

Sample