



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Factor IX

Concentration: 1.44 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.347 ml

Activity: N/A

Total protein: 100.00 units

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: < -60 °C

Expiration Date: 10/2009

PO: N/A

The Bovine Factor IX is >94% pure on SDS-PAGE gels and shows no reduction upon incubation with 2-mercaptoethanol. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,
Quality Assurance

Sample