



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email: info@enzymeresearch.**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Bovine Factor IXa alpha

**Concentration:** 2.0 mg/ml

**Aliquots:** 1 x 100.00 units

**Volume/aliquot:** 0.250 ml

**Activity:** N/A

**Total protein:** 100.00 units

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 05/2006

**PO:** N/A

The Bovine Factor IXa alpha was prepared from homogeneous Bovine Factor IX by activation with Russell's Viper Venom. This RVV-X is removed after activation. All bovine material is purified from U.S. cows < 30 months of age.

Thank You,  
**Quality Assurance**

Sample