



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Factor IXa beta

Concentration: 1.86 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.269 ml

Activity: 3096.00 PEU/mg

Total protein: 100.00 units

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 03/2009

PO: N/A

The Bovine Factor IXa beta was prepared from homogenous Bovine Factor IX by activation with Bovine Factor XIa. This Bovine Factor XIa was removed after activation. Complete activation was observed on 10% SDS-PAGE gels. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,
Quality Assurance

Sample