



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Bovine Factor Xa

**Concentration:** 1.43 mg/ml

**Aliquots:** 1 x 100.00 units

**Volume/aliquot:** 0.699 ml

**Activity:** 221.13 IU/mg

**Total protein:** 100.00 units

**Buffer:** 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 07/2009

**PO:** N/A

The Bovine Factor Xa was prepared from homogeneous Bovine Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation was observed on SDS-PAGE. All bovine material is from USDA inspected herd. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,  
**Quality Assurance**

Sample