



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Factor Xa

Concentration: 1.43 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.699 ml

Activity: 221.13 IU/mg

Total protein: 100.00 units

Buffer: 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 07/2009

PO: N/A

The Bovine Factor Xa was prepared from homogeneous Bovine Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation was observed on SDS-PAGE. All bovine material is from USDA inspected herd. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,
Quality Assurance

Sample