



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Bovine Prothrombin gla domainless

**Concentration:** 1.42 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.704 ml

**Activity:** N/A

**Total protein:** 1.00 mg

**Buffer:** 50 mM Tris-HCl/ 150 mM NaCl/ pH  
7.9

**Storage:** < -60 °C

**Expiration Date:** 07/2009

**PO:** N/A

Bovine prothrombin with the gla domain removed. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,  
**Quality Assurance**

Sample