



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email: info@enzymeresearch.**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Bovine Protein C

**Concentration:** 2.32 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.043 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM  
Benzamidine/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 05/2006

**PO:** N/A

The Bovine Protein C is >95% on SDS-PAGE gels and shows total reduction upon incubation with 2-mercaptoethanol. All bovine material is from U.S. cows < 30 months of age.

Thank You,  
**Quality Assurance**

Sample