



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Activated Protein C

Concentration: 0.92 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 1.087 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 07/2009

PO: N/A

The Activated Bovine Protein C was activated from homogeneous protein C by activation with purified alpha-thrombin. The alpha-thrombin is removed after activation. Complete activation is observed on SDS-PAGE gels. All bovine material is from U.S. cows </+ 36 months of age.

Thank You,
Quality Assurance

Sample