



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Plasminogen

Concentration: 0.96 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 1.042 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 50 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 12/2006

PO: N/A

Bovine Plasminogen is purified from bovine plasma using a lysine-sepharose column. The protein is >94% as judged by SDS PAGE. All bovine material is from U.S. cows < 30 months of age.

Thank You,
Quality Assurance

Sample