



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Thrombin

Concentration: 1.60 mg/ml

Aliquots: 1 x 1000.00 units

Volume/aliquot: 0.294 ml

Activity: 2132.00 NIH u/mg

Total protein: 1000.00 units

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/
0.1% PEG-8000/ pH 6.5

Storage: < -60 °C

Expiration Date: 03/2008

PO: N/A

The Bovine Thrombin was activated from homogeneous Bovine Prothrombin by activation with Factor Va, Factor Xa and phospholipid. These activating enzymes were removed after activation. Complete activation was observed on SDS-PAGE. All bovine material is from U.S. cows < 30 months.

NOTE: Activity is Specific Activity

Thank You,
Quality Assurance

Sample