



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Canine Fibrinogen

Concentration: 23.77 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.042 ml

Activity: 100.00 % Clottable

Total protein: 1.00 mg

Buffer: 20 mM sodium citrate-HCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 05/2007

PO: N/A

The Canine Fibrinogen was purified from fresh frozen canine plasma using a combination of salt precipitations and column chromatography. The purified protein is homogeneous on SDS-PAGE and shows greater than 90% clottability in a thrombin based assay.

Thank You,
Quality Assurance

Sample