



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Canine Factor Xa

Concentration: 2.49 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.040 ml

Activity: 203.00 IU/mg

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 05/2006

PO: N/A

Canine Factor Xa is prepared from homogeneous Canine Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation is observed on SDS-PAGE.

Thank You,
Quality Assurance

Sample