



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human GPIIbIIIa

Concentration: 1.15 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.870 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 0.1%
Triton X-100/ 1 mM CaCl₂/ 0.05%
NaN₃/ 50% glycerol/ pH 7.4

Storage: -20 °C

Expiration Date: 05/2009

PO: N/A

The GPIIbIIIa migrates on gels as follows: GPIIb 136,000 non-reduced and 125,000 reduced. GPIIIa is 97,000 non-reduced and 108,000 reduced. The above protein was purified from Human platelets that were tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample