



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor IXa beta

Concentration: 1.23 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.407 ml

Activity: 10389.00 PEU/mg

Total protein: 100.00 units

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 03/2009

PO: N/A

The Human Factor IXab was activated from homogeneous Human Factor IX by activation with Bovine Factor XIa. This Bovine Factor XIa is removed after activation. Complete activation is observed on SDS-PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample