



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Human Factor VIIa gla domainless

**Concentration:** 0.52 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.192 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 7 mM Tris-HCl/ 52.5 mM sodium  
chloride/ 1.4 mM EDTA/ pH 8.6/ 30%  
glycerol

**Storage:** -20 °C

**Expiration Date:** 07/2008

**PO:** N/A

The Human Factor VIIa gla domainless was prepared from purified Human Factor VIIa using Cathepsin G. The Cathepsin G is removed using affinity chromatography after digestion. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,  
**Quality Assurance**

Sample