



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Human Factor VIIa

**Concentration:** 1.87 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.535 ml

**Activity:** 41870.00 PEU/mg

**Total protein:** 1.00 mg

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 06/2009

**PO:** N/A

The Human Factor VIIa was prepared from homogeneous Human Factor VII using Human Factor XIIa. The Factor XIIa was removed using a Corn Trypsin Inhibitor column. The FVIIa is >95% pure by SDS-PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,  
**Quality Assurance**

Sample