



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HFVIIai

Lot: SAMPLE SAMPLE

Description: Human Factor VIIa Inactivated

Concentration: 1.29 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.078 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 08/2011

PO: N/A

Human Factor VIIa is inactivated using a commercially available chloromethylketone (PPACK II). This inactivator is then removed via dialysis. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample