



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor XIII

Concentration: 4.27 mg/ml

Aliquots: 1 x 0.25 mg

Volume/aliquot: 0.059 ml

Activity: 2342.00 Loewy u/mg

Total protein: 0.25 mg

Buffer: 50 mM Tris/ 100 mM NaCl/ 1 mM
EDTA/ 10 u/mL aprotinin/ 20%
glycerol/ pH 7.5

Storage: < -60 °C

Expiration Date: 07/2008

PO: N/A

The human factor XIII is homogeneous as judged by SDS-PAGE. One plasma equivalent unit (PEU) is equal to ~108 Loewy units. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample