



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor XIIIa

Concentration: 1.27 mg/ml

Aliquots: 1 x 0.25 mg

Volume/aliquot: 0.197 ml

Activity: 1968.00 Loewy u/mg

Total protein: 0.25 mg

Buffer: 50 mM Tris/ 100 mM NaCl/ 1 mM
EDTA/ 10 u/mL aprotinin / 50 uM
DTT/ 50% glycerol/ pH 7.5

Storage: -20 °C

Expiration Date: 07/2008

PO: N/A

Human factor XIII is cleaved with human alpha thrombin. The thrombin is subsequently removed via chromatography. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample