



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor XIa

Concentration: 1.21 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.083 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M
NaCl/ pH 5.3

Storage: < -60 °C

Expiration Date: 06/2009

PO: N/A

The Human Factor XIa was activated from homogeneous Human Factor XI using Human Factor XIIa. This XIIa was removed using a Corn Trypsin Inhibitor column. The protein is >95% pure via SDS-PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample