



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor X

Concentration: 1.33 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.602 ml

Activity: 98.40 PEU/mg

Total protein: 100.00 units

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: < -60 °C

Expiration Date: 03/2009

PO: N/A

The Human Factor X is >95% pure via SDS-PAGE and shows total reduction upon incubation with 2-mercaptoethanol. Human Factor X molecular weight is 58,800 daltons. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample