



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor Xa

Concentration: 1.37 mg/ml

Aliquots: 1 x 100.00 units

Volume/aliquot: 0.584 ml

Activity: 176.00 IU/mg

Total protein: 100.00 units

Buffer: 20mM Bis Tris/ 0.7 M NaCl/ pH 6.0

Storage: < -60 °C

Expiration Date: 06/2009

PO: N/A

The Human Factor Xa was activated from Human Factor X by activation with Russells' Viper Venom. This RVV-X was removed after activation. Complete activation is observed on SDS-PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample