



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Factor Xab gla domainless

Concentration: 1.62 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.617 ml

Activity: 153.00 IU/mg

Total protein: 1.00 mg

Buffer: 20mM Bis-Tris/ 0.1M NaCl / pH 6.0

Storage: < -60 °C

Expiration Date: 06/2009

PO: N/A

The Human Factor Xa beta gla domainless was prepared from HFX-GD using purified RVV-X. This RVV-X was removed after activation by gel filtration. Complete activation was observed on SDS PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample