



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human gamma Thrombin

Concentration: 2.59 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.386 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 10 mM sodium acetate/ 0.7 M
sodium chloride/ pH 6.0-HCl

Storage: < -60 °C

Expiration Date: 01/2009

PO: N/A

The human Gamma-thrombin is >95% pure via SDS-PAGE and shows no reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA.. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample