



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Prekallikrein

Concentration: 1.75 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.571 ml

Activity: 28.00 PEU/mg

Total protein: 1.00 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M
NaCl/ pH 5.3

Storage: < -60 °C

Expiration Date: 07/2008

PO: N/A

The Human Prekallikrein is >95% pure as judged by SDS-PAGE gels and shows no reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample