



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human alpha Thrombin Inactivated

Concentration: 2.92 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.342 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/
0.1% PEG-8000/ pH 6.5

Storage: < -60 °C

Expiration Date: 07/2008

PO: N/A

Human alpha Thrombin treated with an irreversible inhibitor to block the active site. This inhibitor is removed via dialysis. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample