



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Lys-Plasminogen

Concentration: 2.49 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.402 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 50 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 05/2009

PO: N/A

The Human Lys-Plasminogen was prepared from homogeneous glu-plasminogen by batching with a plasmin-sepharose column. This resin was removed after activation. Complete activation was observed by SDS-PAGE. No plasmin activity was detected using the chromogenic assay S-2251. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample