



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Murine Fibrinogen

Concentration: 9.30 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.108 ml

Activity: 91.00 % Clottable

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.15 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 05/2006

PO: N/A

The Murine Fibrinogen was purified from fresh frozen murine plasma using a combination of salt precipitations and column chromatography. The purified protein is homogeneous on SDS-PAGE gels and shows greater than 90% clottability.

Thank You,
Quality Assurance

Sample