



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email: info@enzymeresearch.**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Mouse Thrombin

**Concentration:** 1.09 mg/ml

**Aliquots:** 1 x 0.05 mg

**Volume/aliquot:** 0.046 ml

**Activity:** 2850.00 NIH u/mg

**Total protein:** 0.05 mg

**Buffer:** 50 mM sodium citrate/ 0.2 M NaCl/  
0.1% PEG-8000/ pH 6.5

**Storage:** < -60 °C

**Expiration Date:** 05/2006

**PO:** N/A

The Mouse Thrombin was activated from homogeneous Mouse Prothrombin by activation with Saw Scaled Viper Venom. The venom was removed after activation. Complete activation was observed on SDS-PAGE

Thank You,  
**Quality Assurance**

Sample