



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Factor VIII

Concentration: 1.0 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 10/2009

PO: N/A

Immunogen: Human Factor VIII Clone ID#:012
Purified by Protein G affinity chromatography from cell supernatant.
Antibody binds human factor VIII and human factor VIIIa in solid phase ELISA. Antibody binds human factor VIII A2 domain in immunoblotting techniques.
We suggest hydrating the antibody to original concentration of 1 mg/mL.

Thank You,
Quality Assurance

Sample