



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal to Human Factor X

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 06/2009

PO: N/A

Immunogen: Human Factor X

Antibody binds human factor X in solid-phase ELISA and immunoblotting techniques.

Approximate working concentration of 1-5 ug/mL. Researcher should titer antibody in specific assay.

Reconstitute to original volume with water.

Thank You,
Quality Assurance

Sample