



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Thrombin

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 01/2010

PO: N/A

Antibody binds human thrombin bound to polystyrene and in immunoblotting. Antibody does not bind prothrombin, prethrombin 1 or prethrombin 2, or bovine thrombin, bound to polystyrene.

Recommended use: 1-5 ug/mL. Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory. For research use only.

Reconstitution/Storage: Store lyophilized antibody at -20C. Reconstitute with distilled water. Centrifuge to remove insoluble material. Store reconstituted antibody at 4C or aliquot and store at -20C. Avoid repeated freeze-thaw cycles.

Thank You,
Quality Assurance

Sample