



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Tissue Plasminogen Activator

Concentration: 1 mg/ml

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 09/2008

PO: N/A

Antibody purified by Protein G affinity chromatography from cell supernatant. Antibody binds human tissue plasminogen activator in solid-phase ELISA and immunoblotting techniques. Reconstitute with water to original volume.

Thank You,
Quality Assurance

Sample