



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Monoclonal antibody to Human Vitronectin

**Concentration:** 1.0 mg/ml (before lyophilization)

**Aliquots:** 1 x 0.50 mg

**Volume/aliquot:** 0.500 ml

**Activity:** N/A

**Total protein:** 0.50 mg

**Buffer:** 10 mM Sodium Phosphate/ 0.15 M  
NaCl/ 0.01% Sodium Azide/ pH 7.4

**Storage:** -20 °C

**Expiration Date:** 06/2009

**PO:** N/A

Antibody purified by Protein G affinity chromatography from cell supernatant. Purity determined by SDS-PAGE to be >98%.  
Antibody binds human vitronectin in solid phase ELISA and immunoblotting procedures. Binds plasma vitronectin when coupled to sepharose. Recommended use: 1-5 ug/mL.  
Reconstitute with water to original concentration.

Thank You,  
**Quality Assurance**

Sample