



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Porcine Factor Xa

Concentration: 1.60 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.062 ml

Activity: 201.70 IU/mg

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 01/2007

PO: N/A

The Porcine Factor Xa was activated from homogeneous Porcine Factor X by activation with Russells' Viper Venom. This RVV-X was removed after activation. Complete activation is observed on SDS-PAGE gels.

Thank You,
Quality Assurance

Sample