



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Porcine Fibrinogen

Concentration: 9.04 mg/ml

Aliquots: 1 x 1.00 g

Volume/aliquot: 110.619 ml

Activity: 80.0 % Clottable

Total protein: 1.00 g

Buffer: 20 mM sodium citrate-HCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 07/2007

PO: N/A

The Porcine Fibrinogen was purified from fresh frozen porcine plasma using a combination of salt precipitations and column chromatography. The purified protein is homogeneous on SDS-PAGE and shows greater than 75% clottability in a thrombin based assay.

Thank You,
Quality Assurance

Sample