



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Russells' Viper Venom

Concentration: 2.39 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.418 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: < -60 °C

Expiration Date: 03/2007

PO: N/A

The Russells' Viper Venom is homogeneous as judged by SDS-PAGE and has a high specific clotting activity.

Thank You,
Quality Assurance

Sample