



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Rabbit Fibrinogen

**Concentration:** 14.83 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.067 ml

**Activity:** 100.00 % Clottable

**Total protein:** 1.00 mg

**Buffer:** 20 mM sodium citrate-HCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 03/2007

**PO:** N/A

The Rabbit Fibrinogen was purified from fresh frozen rabbit plasma using a combination of salt precipitations and column chromatography. The purified protein is >94% pure on SDS-PAGE gels and shows greater than 90% clottability.

Thank You,  
**Quality Assurance**

Sample