



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Rabbit Thrombin

**Concentration:** 1.39 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.072 ml

**Activity:** 1324.00 NIH u/mg

**Total protein:** 0.10 mg

**Buffer:** 50 mM sodium citrate/ 0.2 M NaCl/  
0.1% PEG-8000/ pH 6.5

**Storage:** < -60 °C

**Expiration Date:** 12/2006

**PO:** N/A

The Rabbit thrombin is purified from rabbit plasma. Purity is judged by SDS PAGE to be >90%.

Thank You,  
**Quality Assurance**

Sample