



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Rabbit Factor Xa

Concentration: 3.95 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.025 ml

Activity: 148.00 IU/mg

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 12/2006

PO: N/A

Rabbit Factor Xa is prepared from homogeneous Rabbit Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation is observed via SDS-PAGE.

Thank You,
Quality Assurance

Sample