



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Rabbit Factor X

Concentration: 0.9 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.111 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: < -60 °C

Expiration Date: 05/2006

PO: N/A

Rabbit Factor X is purified via salt precipitations along with ion and affinity chromatography. It is judged >94% pure on SDS PAGE.

Thank You,
Quality Assurance

Sample