



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Human TAFI (Procarboxypeptidase-B)

**Concentration:** 0.124 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.806 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 10 mM Hepes/ 0.15 M NaCl/ pH 7.4

**Storage:** < -60 °C

**Expiration Date:** 02/2010

**PO:** N/A

TAFI is a single chain glycoprotein with an apparent mass of 58,000 daltons as determined by SDS-PAGE under reducing conditions. TAFI is purified from normal plasma using immuno-affinity and ion-exchange chromatography.

The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

For Research Purposes Only.

Thank You,  
**Quality Assurance**

Sample